

# Biogel® Tech gloves

Biogel® Tech gloves are designed for use in critical environments. The double-bagged and heat-sealed packaging complies with hygiene protocols and ensures simplicity of use.

## How to unpack your Biogel Tech gloves while ensuring cleanliness

1



Use a plastic apron to protect your clothing from dirt on the cardboard box. This dirt can otherwise be carried onto the plastic bags inside when opening and unpacking the box.



Each transport box contains 200 pairs of gloves. Each of the four compartments in the TRP has 2x25 pairs of gloves packed in heat-sealed double plastic bags\*.

2



Before you unload the plastic bags, make sure your hands are disinfected. Be careful not to touch the plastic bags with contaminated hands.



Place the bags on a clean and disinfected surface.

\*160 pairs for size 9.

3




To remove the outer plastic bag, we recommend using a pair of scissors. Start with the edges to minimize the risk of contamination. Make sure to then fold the plastic outwards without touching the inner plastic bag.

4



Disinfect your hands before you place the 25 pairs (still in the plastic bag) on another clean disinfected surface.



The gloves can be stored in the inner plastic bag or unpacked and placed in a storage box. Put on proper protective clothing and sterile gloves before unpacking the sterile barriers from the inner plastic bag. Again, a pair of scissors are useful to cut the plastic bag at the bottom to allow easy loading of the storage box.



The gloves can now be stored closer to where they will be used.

### Biogel Tech range – Unwrap the next level safety

Each pair of Biogel Tech gloves is wrapped in a high-quality polyethylene inner wrap to keep particle count to a minimum. The wrapping is constructed to stay open when unfolded to ensure aseptic donning of the gloves. The pair is then packed into a sterile barrier made of laminate polyester and low-density polyethylene.

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00  
The Mölnlycke and Biogel trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care group of companies. ©2022 Mölnlycke Health Care AB. All rights reserved. HQIM003293



**Biogel® Tech**  
Precision at every touch

**Mölnlycke®**