

# Treatment of partial thickness burns: A prospective, randomized controlled trial comparing four routinely used burns dressings in an ambulatory care setting

Aggarwala S. *et al.* Journal of Burn Care & Research. 2021;30;42(5):934-943

## Aims

To compare the effectiveness of four dressings for the treatment of adult partial thickness burns, focusing on re-epithelialisation time and cost-effectiveness.

## Method

Prospective, randomised controlled trial  
Adult partial thickness burns patients who met the inclusion criteria were randomised to one of 4 intervention groups:

1. Mepilex® Ag;
2. Biobrane®;
3. Acticoat®;
4. Aquacel® Ag

## Results

No infections were recorded for the course of the study in any treatment arms.

### Healing time

When adjusted for sex, age, smoking status, burn mechanism, TBSA, and first aid adequacy:

- In the Biobrane® group, there was a 26% increase in days to re-epithelialisation compared with Mepilex® Ag (IRR: 1.26, 95% CI: 1.07–1.48,  $P < 0.01$ )
- In the Acticoat® and Aquacel® Ag groups, there was no statistically significant difference in days to re-epithelialisation compared with Mepilex® Ag but a trend in favor to Mepilex® Ag (IRR: 1.12, 95% CI: 0.9–1.3,  $P = 0.24$  and IRR: 1.12, 95% CI: 0.9–1.3,  $P = 0.23$  respectively)

### Cost-effectiveness

Probabilities that Mepilex® Ag was superior (less expensive and more effective) to the other dressings tested:

- **99%** for Mepilex® Ag vs Biobrane®
- **71%** for Mepilex® Ag vs Acticoat®
- **53%** for Mepilex® Ag vs Aquacel® Ag

**Mepilex® Ag was found to be cost-effective in the treatment of partial thickness burns due to a faster rate of re-epithelialisation and a reduction in the cost of dressings compared to Biobrane®, Acticoat® and Aquacel® Ag.**

# To know more about the study

## Outcomes measured

### Primary outcome

Time to wound healing: days to re-epithelialisation

### Secondary outcomes

- Number of outpatient clinic visits
- Pain: assessed during the initial dressing application and each subsequent dressing change using a numeric scale ranging from 0 (no pain) to 10 (extreme pain)
- Nursing experience: assessed by scoring ease of use, ease of application, and ease of removal of dressings. Parameters were measured on the 5-point Likert scale ranging from 1 (very easy) to 5 (very difficult).
- Scar quality: telephone follow-up call at 3 and 6 months after the burn injury (symptoms, subjective outcome parameters)
- Cost-effectiveness of dressings using ICER method:

$$ICER = \frac{(Cost\ Product\ 1\ group) - (Cost\ Product\ 2\ group)}{(Effect\ Product\ 1\ group) - (Effect\ Product\ 2\ group)}$$

## Additional results

131 partial thickness burn wounds in 119 patients were randomised:

- Mepilex® Ag (n = 35);
- Biobrane® (n = 32);
- Acticoat® (n = 37);
- Aquacel® Ag (n = 27).

### Healing time

Dressing	Median days to re-epithelialisation	Interquartile range	Mean (± standard deviation)	P
Biobrane®	11	8.5 – 13	10.8 ± 2.4	0.06
Mepilex® Ag	9	8 – 10	8.9 ± 2.4	
Acticoat®	9	8 – 11	9.6 ± 3.3	
Aquacel® Ag	9	8 – 12	9.6 ± 3.2	

### Cost-effectiveness

Mean consumable cost-saving per patient using Mepilex® Ag vs other dressings			
Comparison*	Mean cost saving (\$AUD)	95% Confidence Interval	P
Mepilex® Ag vs Acticoat®	136	43 – 228	<0.01
Mepilex® Ag vs Aquacel® Ag	2.60	-19 – 24	0.81
Mepilex® Ag vs Biobrane®	148	64 – 233	<0.01

\*Adjusted for total body surface area (%) and mechanism of burn injury.