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Din guide til en mere tryk praksis

Om forfatteren

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Faglige tilknytninger

- Royal College of Obstetricians & Gynaecologists – medlem/akkrediteret underviser
- Medlem af Academy of Medical Educators
- Medlem af OrcID <https://orcid.org/0000-0002-5406-3862>
- Medlem af I.C.S.P (International Council for Surgical Plume), UK Council Member
- Medlem af NAMDET (National Association of Medical Device Educator and Trainers)
- Medlem af Ulster Society of Obstetricians & Gynaecologists (æresmedlem)
- Komitémedlem hos British Standards Institute (ISO)

Publikationer

- Veck. S Eliminating the Hazard of Surgical Plume
Clinical Services Journal March 2021
- Veck. S Smoke Plume – The Risks. National Association of Medical Device Educators & Trainers Journal p21-p22 May 2018
- Veck. S Dorman.G BSCCP Poster May 2016 – Insulated Speculum Use Within A Colposcopy Setting. Are Insulated Vaginal Speculum safe?
- Veck. S Farquharson. R.G
M.R.C.O.G Vignettes in Gynaecology MRCOG Part 1 Supplement
- Veck. S An Introduction to the Principles and Safety of Electrosurgery.Br J Hosp Med 1996. Jan 17–Feb 6; 55 (1-2): 18 19

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Vigtige publicerede artikler

<https://www.hse.gov.uk/research/rrhtm/rr922.htm>

RR922 - Evidence for exposure and harmful effects of diathermy plumes (surgical smoke) – Evidence based literature review

The methods used to dissect tissue and stem blood flow during surgery have changed as technology has developed. Lasers and electro-surgery have become commonplace, so that medical staff in the operating theatre are (potentially) increasingly exposed to the thermal decomposition products of tissues. Variations in ventilation systems and the presence or absence of local exhaust ventilation are likely to influence the extent to which this occurs. A systematic review was carried out to identify existing evidence about surgical smoke (known as diathermy plume) and the potential harm to health care workers exposed in operating theatres. Limited published data were identified, but indicated that dedicated smoke evacuation/extraction devices are effective at reducing the levels of surgical smoke during various surgical procedures, and that correct (close) positioning of smoke evacuation devices to source emissions is likely to be important to the efficiency of surgical smoke removal. The data were insufficient to allow conclusions to be drawn on reported respiratory ill health symptoms linked with surgical smoke exposure.

[https://www.mercyhospital.org.nz/assets/Policies/](https://www.mercyhospital.org.nz/assets/Policies/ElectrosurgicalSmokeEvacuation.pdf)

ElectrosurgicalSmokeEvacuation.pdf

Surgical smoke generated during surgical cases is potentially hazardous and must be captured and filtered through the use of smoke evacuators or in-line filters positioned on suction lines. Surgical smoke (plume) can contain toxic gases and vapours such as benzene, hydrogen cyanide, and formaldehyde along with bio aerosols, dead and live cellular material (including blood fragments), and viruses. At high concentrations, surgical smoke can cause ocular and upper respiratory tract irritation in healthcare workers and can create obstructive visual problems for the surgeon. Surgical smoke has unpleasant odours and has been shown to have mutagenic potential.

www.clinicalservicesjournal.com

Surgical Staff Safety: Going Up in smoke. July 2020

A reader survey has shown that over two-thirds of respondents working in operating theatres are concerned about the effects of surgical smoke on their health, yet only 21% said that their theatres 'always' used smoke evacuation devices when performing electrosurgery or laser treatments. Should their use now become mandatory? Louise Frampton reports.

Journal of Cancer 2019; 10(12):2788-2799

Awareness of surgical smoke hazards and enhancement of surgical smoke prevention among the gynecologists

Yi Liu, Yizuo Song, Xiaoli Hu, Linzhi Yan, and Xueqiong Zhu

Author information Article notes Copyright and License information Disclaimer

Resumé

Surgical smoke is the gaseous by-product produced by heat generating devices in various surgical operations including laser coagulation and loop electrosurgical procedures that often are performed by gynecologists. Surgical smoke contains chemicals, blood and tissue particles, bacteria, and viruses, which has been shown to exhibit potential risks for surgeons, nurses, anesthesiologists, and technicians in the operation room due to long term exposure of smoke. In this review, we describe the detailed information of the components of surgical smoke. Moreover, we highlight the effects of surgical smoke on carcinogenesis, mutagenesis, and infection in gynecologists. Derudover drøftede vi, hvordan man forhindrer kirurgisk røg ved brug af masker med høj filtrering og røgugsystemer samt juridiske retningslinjer for beskyttelsestiltag blandt gynækologerne.

Keywords: Cervical cancer, Cervical intraepithelial neoplasia, Electrosurgery, Smoke, Gynecologist.

Journal of Aerosol Science. 142 (2020) 105512

Morphological Characterization of Particles Emitted from Monopolar Electro Surgical Pencils.

Monopolar electrosurgical pencils are used extensively in surgical operations. With such pencils, electric current passes to the tissue, and as such, electrosurgical pencil operation generates a significant amount of thermal energy, which in turn leads to the generation of electrosurgical smoke (ES). The health risks of ES are dependent on the size distributions as well as the morphologies of the produced particles. To better characterize such particles, in this study we utilized (1) differential mobility analysis with a condensation particle counter (DMA-CPC), (2) an aerodynamic particle spectrometer (APS), (3) DMA-transmission electron microscopy analysis (DMA-TEM), and (4) DMA-aerosol particle mass analysis (DMA-APM) to examine the size distribution and morphologies of particles produced during simulated operation of an electrosurgical pencil (Neptune E-SEP, Stryker Corporation) on bovine, porcine, and ovine tissue. We find that under a variety of operating conditions, ES particles are broadly distributed, with a mode mobility diameter in the 150–200 nm size range, and concentrations well above background levels in the 50nm–5µm size range. We also find that the 'cut' mode of monopolar electrosurgical pencil operation generates higher particle concentrations than the 'coagulate' mode, and that increasing the maximum applied power from 20W to 50W also increases ES particle concentrations. TEM images of mobility selected particles reveal both spherical particles and fractal-like agglomerates in ES; these different particle types are produced under the same operation conditions leading to an externally-mixed, morphologically-complex aerosol. Quantitative analysis of the agglomerate images revealed that agglomerates have an average fractal dimension near 1.93 and that they are structurally similar to agglomerates expected from a diffusion limited cluster aggregation growth mechanism. Despite the presence of both spheres and agglomerates, DMA-APM analysis reveals that all particles have effective densities in the 1000–2000kg m⁻³ range, suggesting that they likely contain inorganic components. Finally, we determined that the collection efficiency of the ES capture suction unit attached to the electrosurgical pencil was >95% for particles in the 50–400nm mobility diameter range.

British Journal of Surgery. BJS May 2020;107:1406-1413

Safe management of surgical smoke in the age of COVID-19

Background: The COVID-19 global pandemic has resulted in a plethora of guidance and opinion from surgical societies. A controversial area concerns the safety of surgically created smoke and the perceived potential higher risk in laparoscopic surgery. Methods: The limited published evidence was analysed in combination with expert opinion. A review was undertaken of the novel coronavirus with regards to its hazards within surgical smoke and the procedures that could mitigate the potential risks to healthcare staff. Results: Using existing knowledge of surgical smoke, a theoretical risk of virus transmission exists. Best practice should consider the operating room set-up, patient movement and operating theatre equipment when producing a COVID-19 operating protocol. The choice of energy device can affect the smoke produced, and surgeons should manage the pneumoperitoneum meticulously during laparoscopic surgery. Devices to remove surgical smoke, including extractors, filters and non-filter devices, are discussed in detail. Conclusion: There is not enough evidence to quantify the risks of COVID-19 transmission in surgical smoke. However, steps can be undertaken to manage the potential hazards. The advantages of minimally invasive surgery may not need to be sacrificed in the current crisis.

Referencer

1. Bree K., et al. (2017). The Dangers of Electrosurgical Smoke to Operating Room Personnel. A Review. *Workplace Health & Safety*, vol. 65, no. 11.
2. <https://www.cdc.gov/niosh/topics/healthcarehps/smoke.html>.
3. Andréasson S.N., et al. (2009). Peritonectomy with high voltage electrocautery generates higher levels of ultrafine smoke particles. *Eur J Surg Oncol*.Jul;35(7):780-4.
4. Rioux M. et al. (2013). HPV positive tonsillar cancer in two laser surgeons: case reports. *J Otolaryngol Head Neck Surg*. 2013;42:54.
5. Alleviating the dangers of surgical smoke. *Quick Safety*, Dec 2020, Issue 56.

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