


Mextra[®] Superabsorbent

Superabsorbent dressing with a fluid-repellent backing film


- High fluid retention capacity¹
- Increased fluid absorption capacity²
- Soft and conformable¹
- Protects the peri-wound area¹




How to use Mextra® Superabsorbent



1. Clean the wound area and select a suitable size of Mextra Superabsorbent.



2. Apply Mextra Superabsorbent directly on the wound area with the white side of the dressing onto the wound. For best results, the pad should overlap the dry surrounding skin by at least 2 cm. Mextra Superabsorbent can also be used together with a wound contact layer, such as Mepitel® and Mepilex® Transfer.



3. Secure Mextra Superabsorbent with a suitable bandage such as Tubifast® bandage or a fixation tape. It can also be used under a compression bandage where appropriate.

Key benefits of Mextra Superabsorbent

- Minimizes the risk of leakage and maceration¹
- Minimizes the need for frequent dressing changes¹
- Minimizes treatment costs¹
- Maintains structure and integrity without becoming bulky upon exudate absorption¹
- Protects against fluid strike-through¹
- Promotes patient comfort¹
- Has protease modulating properties^{1,3} and provides a conducive environment for healing to proceed

Areas of use

Mextra Superabsorbent is intended for use on moderately to highly exuding wounds. Mextra Superabsorbent absorbs wound exudate through the wound contact layer, retains the exudate within the core and thereby reduces the risk of maceration. The fluid-repellent nonwoven backing acts as an exudate barrier and prevents exudate strike-through.

Frequency of change

The wound and Mextra Superabsorbent should be inspected regularly. The dressing should be changed based on the condition of the wound or before the dressing has become saturated. Mextra Superabsorbent may be left on the wound for a maximum of seven days.



Mextra Superabsorbent assortment (Sterile single packed)

Art. no	Size cm	Size inch	Pcs/shelf cont.	Pcs/transp. cont.
610000	12,5 x 12,5	5 x 5	10	50
610100	12,5 x 17,5	5 x 7	10	40
610200	12,5 x 22,5	5 x 9	10	90
610300	17,5 x 22,5	7 x 9	10	40
610400	22,5 x 27,5	9 x 11	10	40
610500	22,5 x 32,5	9 x 13	10	50
610600	22,5 x 42,5	9 x 17	10	50

NOTE: Not all articles are available in every country. Please contact your local Mölnlycke Health Care representative for information about articles available in your country.

REFERENCES: 1. Tickle, J., Fletcher, J. Mextra Superabsorbent made easy. Wounds UK 2012; 8(4): 1-4. 2. Mölnlycke Health Care Laboratory report no: PD-522474, 2016 (unpublished). 3. Wiegand, C., Hipler, U.C. In vitro evaluation of MMP binding capacity of a superabsorbent dressing and the reduction of collagenase activity. Poster presentation at European Wound Management Association conference, Copenhagen, Denmark, 2013. 4. Eming, S., Smola, H., Hartmann, B., Malchau, G., Wegner, R., Krieg, T., Smola-Hess, S. The inhibition of matrix metalloproteinase activity in chronic wounds by a polyacrylate superabsorber. Biomaterials 2008; 29(19): 2932-2940

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3c, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00. The Mölnlycke Health Care, Mextra®, Mepilex®, Mepitel®, and Tubifast® trademarks, names and respective logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. © 2016 Mölnlycke Health Care AB. All rights reserved. HQWC0125 -2016/12

Mölnlycke®