

PRODUCT DATA SHEET

MELGISORB PLUS®

Product Name	Melgisorb Plus®					
Product Description	Absorbent alginate dressing					
Indications	<p>Melgisorb Plus is intended for a wide range of moderately to heavily exuding wounds, both infected and noninfected, such as:</p> <ul style="list-style-type: none"> • Pressure sores • Venous and arterial ulcers • Diabetic ulcers • Donor sites • Post operative wounds • Dermal lesions and other external wounds inflicted by trauma. <p>Melgisorb Plus is not indicated for dry wounds, third-degree burns or surgical implantation. If infection is suspected, follow local routines given by health care professionals.</p>					
Product Codes	Art No	Description	Size (Cm)	Pack Count Con/Ret/Trp	Volume Dm3/Trp	Comments
	252000-00	Melgisorb Plus	5x5	1/10/100		
	252200-00	Melgisorb Plus	10x10	1/10/100		
	252240-00	Melgisorb Plus	10x10	1/3/30		ES
	252300-00	Melgisorb Plus	10x12	1/10/100		
	252500-00	Melgisorb Plus	10x20	1/10/101		
	252540-00	Melgisorb Plus	10x20	1/3/30		ES
	252380-00	Melgisorb Plus	10x12	1/16/160		FR
	252580-01	Melgisorb Plus	10x20	1/16/160		FR
	253500-00	Melgisorb Plus, Cavity	3x45	1/5/50		
	253580-01	Melgisorb Plus, Cavity	3x45	1/16/80		FR
Manufacturer	Mölnlycke Health Care AB Box 13080 Se-40252 Göteborg Sweden					
Country Of Origin	China					
Product Classification According To MDD	Class IIb					

CE-Marked	Yes. Notified Body BSI 0086			
Sterilization Method	Gamma			
Shelf Life	3 years			
Product Material	Melgisorb Plus dressings are composed of Calcium Sodium Alginate fibres. Melgisorb Plus is available as flat nonwoven dressings and flat rope.			
Environmental Classification	ISO 14001	Green		
Other info of special interest	Product and packaging material does not contain phthalates, is free of animal derived material and free of natural rubber latex.			
Physical Characteristics	Property	Value	Reference Method	
	Wet strength	1.59 N/cm	LWI036	
	Absorbency	24,6 g/100 cm ²	EN 13726-1 3.2	
	Retains exudate	Retention capacity 7.23 g/g	LWI161	
	Lateral wicking/minimal spreading	Lateral wicking area 8.64 cm ²	LWI162	
Packaging	<p>The purpose of the dispenser box and transport box is to contain and protect material from dust, moisture and damage during transport, handling and storage.</p> <p>Unit package is a sterile barrier for the product; sterility is guaranteed unless inner package is damage or opened prior to use.</p>			
GS1 CODES	ART NO	CON	RET	TRP
	252000-00	7323190024872	7323190022885	7323190022878
	252200-00	7323190024889	7323190022847	7323190022830
	252240-00	7323190042272	7323190025732	7323190025725
	252300-00	7323190037933	7323190037858	7323190037841
	252500-00	7323190024896	7323190022809	7323190022793
	252540-00	7323190042289	7323190025817	7323190025800
	252380-00	7323190086047	7323190086030	7323190086023
	252580-01	7323190091126	7323190084524	7323190084517
	253500-00	7323190024902	7323190022762	7323190022755
	253580-01	7323190091133	7323190084494	7323190084487