

# Detachable EndoRetrieval Pouch

Truly customised laparoscopic solutions



The Detachable EndoRetrieval Pouch is a unique device that provides the simple, quick and safe capture and retrieval of resected tissues and organs. The detachable facility offers the ultimate convenience by providing the surgeon with the ability to delay the retrieval of the captured tissue. The suture stay design provides rapid retrieval when required, plus the ability to continue to use the port for instrument exchange.

- Unique design ensures simple intuitive manipulation
- Self-opening pouch design allows quick and simple specimen capture
- Enhanced metal ring design ensures simple deployment and a large permanent pouch opening for capture of even difficult specimens
- High-strength, tear-resistant, non-porous material
- Maximum convenience is assured due to the detachable delayed retrieval facility
- Available in a range of sizes to provide solutions for the retrieval of all tissue and organ types
- Available both as sterile single-packed product and for inclusion in Mölnlycke® procedure trays



Detachable EndoRetrieval Pouch

Ref. No.	Description	Size	Pouch size	Volume	Pcs/disp.
899102	Detachable EndoRetrieval Pouch	10mm	9 x 14.5cm	250ml	5
899103	Detachable EndoRetrieval Pouch	10mm	13 x 18cm	600ml	5
899112	Detachable EndoRetrieval Pouch	12mm	19 x 23cm	1150ml	5
899104	Detachable EndoRetrieval Pouch	15mm	19 x 23cm	1150ml	5

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, P.O. Box 13080, Gamlestadvägen 3 C, SE-402 52 Göteborg, Sweden. Phone + 46 31 722 30 00  
 The Mölnlycke trademark, name and logo is registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
 © 2017 Mölnlycke Health Care AB. All rights reserved. HQSU0322